

FASTENING systems & MAGNETIC Keys

Designed for **exclusive** use with SALVAFIX restrain systems

<p>code SVF3500</p>	<p>* PLASTIC Lock system</p> <p>Recomended use in :</p> <p>Hospitals Geriatric institutions Home Care Health shops</p> <p>Medium Performance Button</p> <p>Intended for applied forces limited till 150 KGS</p> <p>.- Top & Bottom Body = 100% injected Plastic .- Pin = 100% sanded stainless steel</p>	
<p>code SVF3550</p>	<p>* PLASTIC Magnetic KEY</p> <p>For use only with PLASTIC Buttons</p>	

<p>code SVF3510</p>	<p>* ALU Lock system</p> <p>Recomended use in :</p> <p>Psychiatrics institutions Urgencies Prisons Police</p> <p>High Performance Button</p> <p>Intended for applied forces limited till 300 KGS</p> <p>.- Bottom Body = 100% anodized reinforced ALUMINIUM .- Top Body = 100% injected Plastic .- Pin = 100% sanded stainless steel</p>	
<p>code SVF3560</p>	<p>* ALU Magnetic KEY</p> <p>Intended for use with ALU Buttons</p> <p>.- Delivered with double Neodimium Magnets .- Reinforced Body = 100 % aluminium</p>	